

Protocols for Ragworms

Duy Nghia Pham, Inna M. Sokolova

2023-01-07

Contents

1	Homogenization	5
1.1	Materials	5
1.2	Methods	6
2	CEA	7
2.1	Organization	7
3	ETS	9
3.1	Materials	9
3.2	Methods	9
4	Proteins	11
4.1	Materials	11
4.2	Methods	11
5	Carbohydrates	13
5.1	Materials	13
5.2	Methods	13
6	Lipids	15
6.1	Materials	15
6.2	Methods	15
7	MGO	17
7.1	Materials	17
7.2	Methods	17
8	MDA	19
8.1	Materials	19
8.2	Methods	19
9	Enzyme assays	21
9.1	Organization	21

10 GR	23
10.1 Materials	23
10.2 Methods	23
11 GST	25
11.1 Materials	25
11.2 Methods	25
12 TAC	27
12.1 Materials	27
12.2 Methods	27
13 Gene expression	29

Chapter 1

Homogenization

1.1 Materials

- 0.1 M Tris:HCl pH 8.5 containing 153 M MgSO₄, 0.2% w/v Triton X-100, 0.1 mM PMSF:
 - 2.5 M HCl: in 250 mL conical flask pre-filled with 100 mL water add 51.754 mL¹ 37% HCl (d = 1.19 g/mL) using glass pipet. Transfer to 250 mL volumetric flask and fill up with water. Transfer to 250 mL bottle. Store at 25°C.
 - In 1 L beaker add 12.114 g² Tris base (M = 121.14) and 0.0377 g MgSO₄*7H₂O (M = 246.48) then fill with 800 mL water. Add 21-23 mL 2.5 M HCl to reach pH 8.5. Transfer to 1 L volumetric flask and fill up with water. Transfer to 1 L bottle then dissolve 1869 L Triton X100 (d = 1.07 mg/ L). Store at 4°C.
 - 200 mM PMSF stock: in 1.5 mL tube add 0.035 g PMSF (M = 174.2) and 1 mL isopropanol, store at -20 °C.
- 1.5 mL flat³ screw cap tubes pre-filled with 500 mg zirconium beads (250 mg 2 mm beads and 250 mg 1 mm beads).
- 2 mL socket screw cap tubes.
- Avery⁴ freeze-resistant round labels, no. 1-15 for each sample.

¹<https://www.sigmaaldrich.com/DE/en/support/calculators-and-apps/molarity-calculator>

or ²<https://www.sigmaaldrich.com/DE/en/support/calculators-and-apps/mass-molarity-calculator>

³compatible with FastPrep

⁴<https://www.avery-zweckform.com/blanco-etiketten/rund-10-mm>

1.2 Methods

- Quickly cut the head (50 mg) and wrap it in the aluminum foil maintained in liquid nitrogen.
- Record mass of the headless worm, transfer to 15 mL pre-cool glass homogenizer.
- Prepare 1:25 tissue homogenate:
 - in 50 mL tube add 2.5 L PMSF stock to every 5 mL cold buffer (2000-fold dilution) right before use. Invert briefly.
 - Add 2500 L buffer containing PMSF to every 100 mg wet tissue.
- Record mass of the head, transfer to the 1.5 mL tube pre-filled with beads (no. 1). Cool down in liquid nitrogen then temporarily store at -20 °C.
- Surround glass homogenizer by ice, rotate at 200 rpm for ~3 min.
- Aliquot 1.8 mL homogenate to 2 mL tubes (no. 2+) on ice then temporarily store at -20 °C.
- Store tubes (no.1 for gene expression and no. 2+ for other biomarkers in separated cryoboxes) at -80 °C for downstream applications.

Chapter 2

CEA

2.1 Organization

- Transfer one homogenate tube (e.g., no. 2) of all samples to new cryoboxes.
- Thaw out homogenate (max. 24 samples each run, 8 samples each sub-run).
- Vortex and transfer 375 μ L homogenate to 2 mL tube, store at -80 °C for later analysis of lipids (read more in 6).
- Vortex and transfer 300 μ L homogenate to 1.5 mL tube, centrifuge 3,000 $\times g$ for 5 min at 4 °C.
- Transfer supernatant to 1.5 mL tube, use for ETS on the same day (read more in 3) and store the leftover at -80 °C for later analysis of carbohydrates and proteins (read more in 4 and 5).

Chapter 3

ETS

Modified after De Coen and Janssen (1997).

3.1 Materials

- Buffer solution: 0.13 M Tris:HCl pH 8.5 containing 0.3% w/v Triton X-100
 - Add 1.574 g Tris base ($M = 121.14$) to conical flask then fill with ~80 mL water. Add ~3.1 mL 2.5 M HCl to reach pH 8.5. Add 280.4 L Triton X100 ($d = 1.07 \text{ mg/L}$). Dissolve, adjust pH if needed and add water to a volume of 100 mL, aided by volumetric flask. Transfer to DURAN bottle, store cold.
- Substrate solution: 1.7 mM NADH and 0.25 mM NADPH Add 0.0241 g NADH ($M = 709.4$) and 0.0042 g NADPH ($M = 833.35$) to 20 mL buffer solution, store dark cold
- 8 mM INT: 0.2023 g INT ($M = 505.7$) in 50 mL water, store dark cold
- For each batch of 24 samples, take out 10 mL buffer in Falcon tube, 1.3 mL Substrate in 1.5 mL tube, 5.3 mL INT in 50 mL glass bottle
- 0.5 M KCN ($M = 65.12$): 32.6 mg KCN in 1mL water, freeze at -20
- 1 mM rotenone ($M = 394.42$): 0.39 mg rotenone in 1ml ethanol absolute, freeze at -20

3.2 Methods

- Rxn in duplicate wells:
 - 25 L supernatant + 75 L buffer + 23 L KCN + 2 L rotenone + 50 L INT for self-control
 - 25 L supernatant + 75 L buffer + 25 L substrate + 50 L INT (Handystep) for rxn

- Measure kinetically at 490 nm in 10 min with 10s interval, at room temperature (25 °C), shake 5s before reading

Chapter 4

Proteins

Modified after Bradford (1976).

4.1 Materials

- 5 mg/mL BSA stock solution: e.g., 50 mg BSA in 10 mL homogenization buffer
- 2 mg/mL BSA stock solution: 0.5 mL 5 mg/mL BSA + 0.75 mL water
- Bradford reagent

4.2 Methods

- dilute supernatant 8 times in 1.5 mL tube: 10 μ L supernatant + 70 μ L water, on ice
- Prepare standard curve using 2 mg/mL BSA stock solution in 1.5 mL tubes on ice

mg/mL standard	μ L stock	μ L water
2	200	0
1	100	100
0.5	50	150
0.25	25	175
0.125	12.5	187.5
0.0625	6.25	193.75
0	0	200

- Add 200 μ L Bradford reagent using handystep to 10 μ L diluted supernatant or standards in wells
- Shaking 300s (5 min) and measure in triplicate at 595 nm

Chapter 5

Carbohydrates

Modified after Masuko et al. (2005).

5.1 Materials

- Concentrated H₂SO₄, store dark room
- 5% phenol: 2 mL 90% phenol + 34 mL water in Falcon tube, vortex then transfer to 50 mL glass bottle, store dark
- 100 mM glucose stock solution: e.g., 180.16 mg in 10 mL water, store cold
- 2 mM glucose stock solution: 200 μ L stock solution + 9800 μ L water, store cold

5.2 Methods

- dilute supernatant 8 times in 2 mL tube: 25 μ L supernatant + 175 μ L water = 200 μ L
- Prepare standard curve using 2 mM stock solution in 2 mL tubes: 2, 1, 0.5, 0.25, 0.125, 0.0625, 0 mM (200 μ L in final)

mM standard	μ L stock	μ L water
2	200	0
1	100	100
0.5	50	150
0.25	25	175
0.125	12.5	187.5
0.0625	6.25	193.75
0	0	200

- Medium vortex samples and standards
- Add 500 μ L concentrated H_2SO_4 and 100 μ L 5% phenol both using Handystep
- VORTEX FOR A WHILE and heat at 90 $^\circ\text{C}$ for 5 min, cool down for 5 min then vortex.
- measure 200 μ L in triplicate at 492 nm. Shake 5s before reading.

Chapter 6

Lipids

Modified after Inouye and Lotufo (2006).

6.1 Materials

- Chloroform: methanol 1:1 v/v: mix 150 mL chloroform and 150 mL methanol in a 500 mL bottle, store dark room, pour to 100 mL bottle to use
- Concentrated H₂SO₄, store dark room
- 5 mg oil/mL acetone stock solution: 5.7 L oil ($d = 0.87 \text{ mg/L}$) in 1 mL acetone¹
- Vanillin reagent: dissolve 750 mg vanillin in 125 mL hot water then add 500 mL 85% H₃PO₄ in a 1 L bottle, store dark room, pour to 250 mL bottle to use

6.2 Methods

- Prepare standard curve using stock solution in 1.5 mL tubes: 5, 2.5, 1.25, 0.625, 0.3125, 0 mg/mL

mg/mL standard	uL stock	uL acetone
5	200	0
2.5	100	100
1.25	50	150
0.625	25	175
0.3125	12.5	187.5

¹Pre-wetting 1 mL tip can prevent acetone dropping

mg/mL standard	uL stock	uL acetone
0	0	200

- Add 375 μ L homogenate to 2 mL tube²
- thawing homogenate, medium vortex, add 1.5 mL Chloroform: methanol, vortex and incubate 10 min
- Centrifuge 3000 g for 5 min at 4 $^{\circ}$ C, wait some time for methanol going up
- transfer 200 μ L of low chloroform phase ³ to 1.5 mL tube
- spin down standards and samples
- evaporate solvent in standards and samples by heating at 90 $^{\circ}$ C for 10 min (open cap)
- Add 100 μ L concentrated H₂SO₄ using Handystep, vortex twice and heat at 100 $^{\circ}$ C for 10 min (close cap)
- After cooling, transfer 50 μ L to 1.5 mL microtubes tube already filled with 1 mL vanillin reagent, vortex
- After 5 min, measure 200 μ L in triplicate at 490 nm

²already done when preparing for ETS assay

³wait for methanol getting out of the tip

Chapter 7

MGO

Modified after Mitchel and Birnboim (1977).

7.1 Materials

- 0.1 M borax: in 50 mL tube add 1.5255 g borax ($M^1 = 381.37$) and 40 mL water. Make another one. Rotate 20 min. Transfer to 100 mL bottle. Store in dark at 25°C².
- 0.2 M Girard's reagent T: in 50 mL tube add 1.3411 g Girard T ($M = 167.64$) and 40 mL water. Rotate 5 min. Transfer to 50 mL bottle. Store in dark at 4°C.
- 1232 uM methylglyoxal stock solution: in 50 mL tube add 10 uL 6.16 M methylglyoxal ($M = 72.06$, 40% w/w in water, $d = 1.11$) and 50 mL homogenization buffer³ (5000-fold dilution). Rotate 5 min. Store at 4°C.

7.2 Methods

- Thaw out homogenate (max. 24 samples each run). **Vortex** and transfer 250 uL to 1.5 mL tube. Centrifuge at $10,000 \times g$ for 15 min at 25°C. Transfer 200 uL supernatant to 1.5 mL tube.
- Prepare 200 uL methylglyoxal standards in 1.5 mL tube by stock dilution.

uM standard	uL stock	uL homogenization buffer
1232	200	0
616	100	100

¹borax decahydrate

²avoid lower temperature to prevent precipitation

³0.1 M Tris:HCl pH 8.5 containing 153 uM MgSO₄, 0.2% w/v Triton X-100, no PMSF and EDTA

uM standard	uL stock	uL homogenization buffer
308	50	150
154	25	175
77	12.5	187.5
38.5	6.25	193.75
0	0	200

- Add 250 uL borax and 50 uL Girard T to samples and standards using repeating pipette.
- **Vortex** and heat at 60°C⁴ for 10 min. Cool to 25°C for 5 min.
- Centrifuge at 10,000 $\times g$ for 5 min at 25°C⁵.
- Measure⁶ absorbance of 200 μ L in **duplicate** at 325 nm. Shake 5s before reading.

⁴trials at 30°C and 90°C showed lower absorbance

⁵may need another centrifuge for standards; small pellet could appear in samples but not in standards

⁶pre-wetting the tip helps minimize Triton bubbles and increase homogeneity in samples (since supernatant is taken directly for measurement without transferring to new tube and vortexing)

Chapter 8

MDA

Modified after Buege and Aust (1978).

8.1 Materials

- 15% w/v TCA - 0.375% w/v TBA - 0.25 M HCl reagent: Add 0.1875 g TBA in 50 mL tube. Add 34 mL 20% w/w TCA (density 1.1 g/mL, equivalent to 22% w/v). Add 5 mL 2.5 M HCl. Add 10 mL water to a volume of nearly 50 mL. Wrap by aluminum, rotate for 1h, store dark room. Rotate for 5 min before use each day.

8.2 Methods

- Thaw out homogenate (max. 22 samples each run). **Vortex** and transfer 200 uL to 1.5 mL tube **with snap cap with hole made by needle**. Use 200 uL homogenization buffer for control in **duplicate**.
- Add 500 uL reagent to samples and controls (24 tubes in total), **NOT vortex**.
- Heat at 90 °C for 20 min. Cool for 15 min.
- Centrifuge at 10,000 $\times g$ for 3 min at 25 °C. Rotate tubes 180° then centrifuge one more time with same settings.
- Allow air to escape (especially in control). Measure absorbance of 200 L in **triplicate** at 530 nm. Shake 5s before reading.

Chapter 9

Enzyme assays

9.1 Organization

- Thaw out homogenate (max. 22 samples each run, 11 samples each sub-run).
- Vortex and transfer 300 μ L homogenate to 1.5 mL tube, centrifuge 6,000 $\times g$ for 10 min at 4 $^{\circ}\text{C}$.
- Transfer supernatant to 1.5 mL tube, use for GR on the same day (read more in 10) and store the leftover at -80 $^{\circ}\text{C}$ for later measurement of GST (read more in 11).

Chapter 10

GR

Modified after Mannervik (1999).

10.1 Materials

- 1 mM NADPH: in 5 mL tube add 4.2 mg NADPH ($M = 833.4$) and 5 mL homogenization buffer. Store at 4°C.
- 10 mM GSSG: in 5 mL tube add 30.6 mg GSSG ($M = 612.6$) and 5 mL water. Store at 4°C.
- 1% BSA: in 5 mL tube add 50 mg BSA and 5 mL homogenization buffer. Store at 4°C.

10.2 Methods

- Master mix for 24 rxn: in 5 mL tube add 3120 uL homogenization buffer, 520 uL NADPH, 520 uL GSSG, and 520 uL BSA (120:20:20:20). **Vortex.**
- Rxn in duplicate wells:
 - 180 uL master mix + 20 uL homogenization buffer for control
 - 180 uL master mix + 20 uL supernatant for sample
- Measure kinetically at 340 nm in 15 min with 30s interval, at room temperature (25 °C), shake 5s before reading. Use 21/31 points for slope.

Chapter 11

GST

Modified after Habig et al. (1974).

11.1 Materials

- 5 mM GSH: in 5 mL tube add 7.7 mg GSH ($M = 307.33$) and 5 mL homogenization buffer. Use fresh only.
- 5 mM CDNB:
 - Prepare 100 mM CDNB stock: e.g., in 15 mL tube add 286.9 mg CDNB ($M = 202.55$) and 14.16 mL absolute ethanol. Store at -20°C .
 - Dilute 20X: in 2mL tube add 50 uL stock and 950 uL absolute ethanol.

11.2 Methods

- Master mix for 24 rxn: in 5 mL tube add 4160 uL homogenization buffer, 260 uL GSH, and 260 uL CDNB (160:10:10). **Vortex**.
- Rxn in duplicate wells:
 - 180 uL master mix + 20 uL homogenization buffer for control
 - 180 uL master mix + 20 uL supernatant for sample
- Measure kinetically at 340 nm in 15 min with 30s interval, at room temperature (25°C), shake 5s before reading. Use 21/31 points for slope.

Chapter 12

TAC

Modified after Re et al. (1999).

12.1 Materials

- 63.63 mM potassium persulfate: in 5 mL tube add 86 mg $K_2S_2O_8$ ($M = 270.32$) and 5 mL water. Vortex.
- 7.29 mM ABTS: in 15 mL tube add 40 mg ABTS salt ($M = 548.68$) and 10 mL water. Vortex.
- ABTS radical cation stock solution: add 0.4 mL $K_2S_2O_8$ to 10 mL ABTS (final conc. 2.45 mM $K_2S_2O_8$ + 7 mM ABTS). Vortex. Store in dark at 25°C for 24h before use.
- 10 mM Trolox: in 5 mL tube add 12.5 mg Trolox ($M = 250.29$) and 5 mL homogenization buffer. Rotate 5 min. Use fresh.

12.2 Methods

- Thaw out supernatant (max. 40 samples each run). **Vortex** before use.
- Prepare 200 uL Trolox standards in 1.5 mL tube by stock dilution. **Vortex** and keep on ice.

uM standard	uL stock	uL homogenization buffer
0	0	200
500	10	190
1000	20	180
1500	30	170
2000	40	160
2500	50	150

uM standard	uL stock	uL homogenization buffer
3000	60	140

- Dilute¹ ABTS^{•+} 10X: in 50 mL tube add, e.g., 2.5 mL ABTS^{•+} stock solution and 22.5 mL water. Vortex.
- Rxn in duplicate wells: 10 uL standards or samples + 200 uL diluted ABTS^{•+} using repeating pipette.
- Measure endpoint at 734 nm at room temperature (25 °C), shake 300s (5 min) before reading.

¹check Abs beforehand with 3000 and 0 uM standards to fit the reader capacity. Adjust by adding stock solution or water

Chapter 13

Gene expression

updating

Bibliography

- Bradford, M. M. (1976). A rapid and sensitive method for the quantitation of microgram quantities of protein utilizing the principle of protein-dye binding. *Analytical Biochemistry*, 72(1-2), 248–254. [https://doi.org/10.1016/0003-2697\(76\)90527-3](https://doi.org/10.1016/0003-2697(76)90527-3)
172736 citations (Crossref) [2022-08-26].
- Buege, J. A., & Aust, S. D. (1978, January 1). [30] Microsomal lipid peroxidation. In S. Fleischer & L. Packer (Eds.), *Methods in Enzymology* (pp. 302–310, Vol. 52). Academic Press. [https://doi.org/10.1016/S0076-6879\(78\)52032-6](https://doi.org/10.1016/S0076-6879(78)52032-6)
- De Coen, W., & Janssen, C. (1997). The use of biomarkers in *Daphnia magna* toxicity testing. IV. Cellular Energy Allocation: A new methodology to assess the energy budget of toxicant-stressed *Daphnia* populations. *Journal of Aquatic Ecosystem Stress and Recovery*, 6(1), 43–55. <https://doi.org/10.1023/A:1008228517955>
311 citations (Crossref) [2022-08-26].
- Habig, W. H., Pabst, M. J., & Jakoby, W. B. (1974). Glutathione S-transferases. The first enzymatic step in mercapturic acid formation. *The Journal of Biological Chemistry*, 249(22), 7130–7139. [https://doi.org/10.1016/S0021-9258\(19\)42083-8](https://doi.org/10.1016/S0021-9258(19)42083-8)
- Inouye, L. S., & Lotufo, G. R. (2006). Comparison of macro-gravimetric and micro-colorimetric lipid determination methods. *Talanta*, 70(3), 584–587. <https://doi.org/10.1016/j.talanta.2006.01.024>
40 citations (Crossref) [2022-08-26].
- Mannervik, B. (1999). Measurement of Glutathione Reductase Activity. *Current Protocols in Toxicology*, 00(1), 7.2.1–7.2.4. <https://doi.org/10.1002/0471140856.tx0702s00>
27 citations (Crossref) [2022-08-26].
- Masuko, T., Minami, A., Iwasaki, N., Majima, T., Nishimura, S.-I., & Lee, Y. C. (2005). Carbohydrate analysis by a phenol–sulfuric acid method in microplate format. *Analytical Biochemistry*, 339(1), 69–72. <https://doi.org/10.1016/j.ab.2004.12.001>
1273 citations (Crossref) [2022-08-26].

- Mitchel, R. E. J., & Birnboim, H. C. (1977). The use of Girard-T reagent in a rapid and sensitive method for measuring glyoxal and certain other -dicarbonyl compounds. *Analytical Biochemistry*, 81(1), 47–56. [https://doi.org/10.1016/0003-2697\(77\)90597-8](https://doi.org/10.1016/0003-2697(77)90597-8)
51 citations (Crossref) [2022-08-26].
- Re, R., Pellegrini, N., Proteggente, A., Pannala, A., Yang, M., & Rice-Evans, C. (1999). Antioxidant activity applying an improved ABTS radical cation decolorization assay. *Free Radical Biology and Medicine*, 26(9), 1231–1237. [https://doi.org/10.1016/S0891-5849\(98\)00315-3](https://doi.org/10.1016/S0891-5849(98)00315-3)